



9781 Green Park Industrial Drive  
P.O. Box 4358  
St. Louis, Missouri 63123

**LOCAL**  
314-533-3440

**TOLL-FREE**  
800-325-8921

**FAX**  
314-533-7331

**EMAIL**  
Chris@OhlendorfApplianceLab.com

## CLEAR SOLUTIONS PRESCRIPTION SHEET

Dr: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ St: \_\_\_\_\_ Zip: \_\_\_\_\_

Phone: \_\_\_\_\_ Dr. Email: \_\_\_\_\_

Pt. Name: \_\_\_\_\_ Date Wanted: \_\_\_\_\_

☐ Upper ☐ Lower ☐ Please Call Dr. to Discuss Case

\_\_\_\_\_ **Clear Hold Retainer** – Invisible/Essix Retainer to hold the teeth in their current position

\_\_\_\_\_ **Clear Hold Package** – 2 Upper & 2 Lower Clear Hold Retainers & Set of 3D printed models

\_\_\_\_\_ **Clear Express** – 1-5 Clear Tray System to align the upper or lower 4 anterior teeth  
Includes IPR and Movement Guide

\_\_\_\_\_ **Clear Force** – Full Arch Movement up to 20 Trays

Alignment Needed:

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8

\_\_\_\_\_ **Special Instructions:**

Signature: \_\_\_\_\_ License #: \_\_\_\_\_

### How to Send Us a Case

1. Please make a copy of this page to use as a prescription sheet and a mailing label.
2. Fill out the prescription sheet
3. Ship us upper and lower models or PVS impressions with the prescription sheet
4. Tape the mailing label on the outside of the box
5. Deliver package to the U.S. Post Office

Additional mailing supplies can be downloaded by visiting our website:  
[www.OhlendorfApplianceLab.com](http://www.OhlendorfApplianceLab.com)

FROM




**BUSINESS REPLY LABEL**  
FIRST CLASS PERMIT NO. 3735 ST. LOUIS, MO.

POSTAGE WILL BE PAID BY ADDRESSEE



P.O. Box 4358  
St. Louis, MO 63123

**NO POSTAGE  
NECESSARY  
IF MAILED  
IN THE  
UNITED STATES**



## OR SEND US DIGITAL SCANS

*We accept scan from  
the following systems*

**Itero** **CareStream**  
**3Shape** **3M Lava**  
**Sirona** **.stl files**

### STEP 1

Send us Upper and Lower Scans

### STEP 2

Write out instructions in the Notes Section

*If you are unsure how to find us in your system, please call us at  
1-800-325-8921 and we will guide you through the process.*

# CLEAR SOLUTIONS™

## INVISIBLE TRAY SYSTEMS FOR

Individual Tooth  
Movement

Full Arch  
Movement

Retention



### CLEAR EXPRESS™

1 - 5 Tray System to Align Four Anterior Teeth

### CLEAR FORCE™

Full Arch Comprehensive Alignment

### CLEAR HOLD™

Invisible/Essix Retainer

### CLEAR HOLD PACKAGE™

Two Sets of Upper and Lower Clear Hold Retainers



**OHLENDORF**  
**APPLIANCE LABORATORY**

*Working With You to Deliver Smiles ... Industry Best Solutions Since 1933*

(800) 325-8921 • [www.OhlendorfApplianceLab.com](http://www.OhlendorfApplianceLab.com)



## Clear Express™

- 1-5 Invisible Tray system
- Aligns the 4 anterior teeth
- Minor IPR needed to create space to align teeth
- No case is too small
- Pay per tray system
- Correct minor relapse and rotations

### What you will receive:

- ✓ Complete set of trays to finish the case
- ✓ IPR Guide
- ✓ Individual Tooth Movement Guide
- ✓ Final 3D Printed Model



## Clear Hold™

- Invisible / Essix type retainer
- Worn on the upper and/or lower arch
- Fabricated on a BioStar machine for superior fit
- High quality material that last longer and resist staining
- Fit over a fixed lingual bonded retainer

## Clear Force™

- Clear Tray System for Full Arch Movement
- Comprehensive system to move all the teeth
- Correct larger amounts of crowding and rotations
- Most cases average 12 to 16 trays per arch

### What you will receive:

- ✓ All the trays to complete the case
- ✓ Attachment Template if needed
- ✓ Treatment Movement Video
- ✓ Smile Summary for IPR and attachment placements
- ✓ "Chewies" for better seating

**Our expert technicians are available to help you determine the best Clear Solution for your next case**

**CALL US AT (800) 325-8921**

## Clear Hold Package™

- 2 sets of Upper Clear Hold Retainers
- 2 sets of Lower Clear Hold Retainers
- A set of Upper and Lower 3D printed models

Order your retainers as a package and **save 50%** off our individual retail price



### OUR GUARANTEE TO YOU:

- ✓ If you need assistance, you can talk directly to one of our experienced technicians at our Lab in St. Louis about your case
- ✓ We will never sell to the public directly
- ✓ No case is too small
- ✓ You will receive the same high quality personal service no matter how big or small your account is with us
- ✓ Our pricing structure makes our services more reasonable than other companies so more of your patients can afford to have the beautiful smile they deserve



### **Not sure where to begin?**

We offer Free Evaluations for both Clear Express and Clear Force cases. Simply send us your case, and one of our expert Technicians will contact you to:

- Let you know how many trays may be needed
- An estimate of the Lab Cost
- Which system would work best for your patient

If you have questions about a case or need assistance, please call our Appliance Hot Line at 1-800-325-8921.

We will answer your questions and provide the excellent Technical Support you are accustomed to at the Ohlendorf Appliance Laboratory.